

Information for Use

Weighted Tactile Cushion



SKU: WCUS

WARNINGS:



Keep away
from fire



Adult
supervision
required



Do not
use while
sleeping



Use one
weighted
product at
a time

Important:

- Recommended age 8+
- Always ensure the user can independently remove the Weighted Tactile Cushion.
- Remove the Weighted Tactile Cushion at the first signs of discomfort or distress.
- Check before use, discontinue using if showing signs of wear and tear.
- The Weighted Tactile Cushion has a washable, removable cover.
- Do not alter, adapt or wash outside of the care instructions.

PRODUCT OUTER	PRODUCT INNER	PRODUCT FILLING
100% Polyester Minky Dimple Fleece	65% Polyester / 35% Cotton	Polyester staple fibre and Polyethylene HDPE beads
<p>Care Instructions – Cover Only</p> <p>Machine washable 30°C. Use a gentle wash and low spin cycle. Do not use fabric softener or bleach. Air dry. Do not tumble dry.</p>		

Good to know:

- Dimensions - (L) 50cm x (W) 32cm.
- Weight - 2kg.

Good for:



Calming and Relaxing.
Reducing the need
for repetitive sensory
seeking



Improving state of
wellbeing.
Enabling self-regulation



Improving concentration,
focus & attention.
Improving body
awareness



Safe and non-invasive
therapeutic solution

How to use:

- As a guide, the Weighted Tactile Cushion should be used for approximately 20-40 minutes at a time and then removed for the same period before they are used again.
- Ensure that the user is consistently comfortable and content with the deep pressure they experience.

Ideal for:



Education or
Home Study



Dining



Chill-Out Time



Home or
Residential Setting



Travelling



Manufactured by Sensory Direct (UK) Limited

Unit 65c, Blackpole Trading Estate West, Worcester, WR3 8TJ, UK

☎ 01905 670 500 ✉ info@sensorydirect.com 🌐 www.sensorydirect.com