

Information for Use

Therapy Balls / Peanut Balls / Sitting Balls

SKU: BRG / TSA / PNA / PNM / SNH

We recommend consultation with an Occupational Therapist or other suitably qualified person prior to use.

WARNINGS:



Adult supervision required



Contains small parts - choking hazard



Do not use weights when using this product



This product is not a toy



Do not exceed the maximum weight limit

Important:

- **Therapy Balls** – Maximum weight limit: 114kg.
- **Peanut Balls** – Maximum weight limit: 110kg.
- **Sitting Balls** – Maximum weight limit: 130kg.
- Check the Ball for damage before each use. If the Ball is damaged in any way do not use.
- If the Ball develops a puncture do not repair and use.

Good to know:

- Made from PVC – Free from BPA, phthalates, latex and rubber.
- Burst resistant – In case of a puncture, the Ball will deflate slowly.

Inflating the Ball:

- Remove the plug from the Ball and insert a pump to inflate (not included). You can use a hand pump, foot pump or electric pump. If necessary, use the special adaptor provided.
- The Ball should be inflated at room temperature.
- DO NOT EXCEED the Ball diameter size when inflating.

Care of the Ball:

- Clean with warm water and neutral soap.
- Do not use brushes, sponges or chemicals to clean the Ball as this may damage the surface.
- Keep the Ball away from sharp objects or surfaces that may damage it.
- Keep the Ball away from direct heat.
- Do not expose the Ball to direct light for long periods of time.

Made in Taiwan