

Portable Massage Cushion

SKU: VSP45BLK

WARNINGS:

- KEEP AWAY FROM FIRE.
- Adult supervision required.
- Do not use during pregnancy.
- Do not use on bruised or broken skin.
- Do not immerse in water.

Important:

- This is not a toy.
- Do not mix alkaline, standard (carbon zinc) or rechargeable batteries.
- Batteries should be replaced by a responsible adult.

How to use:

- Ensure the cushion is turned off before replacing batteries.
- Unzip the back and insert the batteries into the compartment.
- In the event of pain or nausea when using, discontinue use immediately.
- Recommended usage is a maximum of 15 minutes at a time.

Care instructions:

Made from Polycotton and filled with 100% Polyurethane foam.

- Do not machine wash.
- Clean with a lightly damp cloth - outer only.
- Ensure product is dry before storage
- Store in cool, dry environment
- Remove batteries before cleaning.

Technical information:

- Voltage / Power - 4 x AA batteries (not included).
- Dimensions - 33cm x 30cm x 10cm
- Net Weight - 0.22kg.



Made in China

Imported by Bing International Ltd, Unit 14 Elgar Business Centre, WR2 6NJ on behalf of Sensory Direct (UK) Ltd, Unit 65c, Blackpole West, WR3 8TJ, UK. www.sensorydirect.com

The UKCA logo, consisting of the letters "UK" stacked above "CA" in a bold, sans-serif font.